

## M-Bond 610

Two-component system

### Description

M-Bond 610 is a permanent, two-component epoxy-phenolic adhesive that cures in 1 hour at 150 °C (302 °F). It is chemically resistant, provides a very thin glue line and will not outgas under vacuum. Curing must begin within 4 hours of application.

### Mixing Ratio, by VOLUME

1:1 (Resin: Curing Agent)

### Curing Schedule

(glue line temperature)

1 hour at 150 °C (302 °F)

2.5 hours at 120 °C (248 °F)

4 hours at 100 °C (212 °F)

### Operating Temperature Range

Short Term: -452° to +700°F (-269° to +370 °C)

Long Term: -452° to +500°F (-269° to +260 °C)

Transducers: to +450°F (+230 °C)

### Physical Properties (as-mixed)

Pot Life = 6 weeks, at 24 °C (75 °F) or

12 weeks at +40 °F (+5 °C)

### Elongation Capabilities

1% at -452 °F (-269 °C)

3% at +75 °F (+24 °C)

3% at +500 °F (+260 °C)

### Clamping Pressure

10 to 70 psi (70 to 480 kN/m<sup>2</sup>)

30 to 40 psi optimum (200 to 275 kN/m<sup>2</sup>)

### Cure Requirements

Recommended Postcure: 2 hours at 50° to 75°F (30° to 40°C) above maximum operating temperature or cure temperature, whichever is higher

High Precision Transducer Postcure: 2 hours at +400° to +450°F (+205° to +230°C) after wiring

### Applications

M-Bond 610 is commonly used to bond multiple samples together and to adhere thinned samples to grids for TEM or FIB observation.

## WARNING!



Refer to the SDS document for additional safety information.

### Instructions:

**DO NOT** open a resin or curing agent bottle before it has reached room temperature!

- 1) Using the included disposable plastic funnel, empty the contents of the bottle labeled "Curing Agent" into the bottle labeled "Adhesive." Discard the funnel when finished.
- 2) Use the included brush applicator cap to securely close the "Adhesive" bottle and thoroughly mix the contents by vigorously shaking it for 10 seconds.
- 3) Mark the bottle in the space provided on the label with the mixing date.
- 4) Allow the freshly mixed adhesive to stand for at least 1 hour before using.

**Storage Note:** Keep the containers closed when not in use, and do not store above 48 °C (120 °F). The shelf life is 9 months from date of manufacture when stored at room temperature (+75°F, +24°C) or 15 months from date of manufacture when refrigerated (+40°F, +5°C). Refrigerated storage is recommended.