



TechCut 5™
PRECISION HIGH SPEED SAW

*A VERSATILE, PROGRAMMABLE
SECTIONING MACHINE DESIGNED TO
CUT A WIDE VARIETY AND SIZE OF MATERIALS*

Quality Products for Metallographic Sample Preparation & Analysis

The **TechCut 5™** precision high speed saw is a versatile, programmable machine designed to cut a wide variety and size of materials.

It automatically sections materials at high speeds, increasing sample throughput.

The microprocessor-based system controls sample feed rate, distance and force, and automatically adjusts feed rate as the cutting condition changes due to varying thickness and/or material differences in the sample.

When sectioning is complete, the table automatically retracts the sample to the home position and stops blade rotation and coolant.

The unique fixturing system allows for easy changes between the T-slot table and the X-axis tables. Both table types offer a variety of convenient table-specific fixture options.

Tables, accessories and consumables are sold separately.



Unparalleled Fixturing Versatility

Item	Description
5-5700-ER	TechCut 5™ with External Reservoir, 100-240 V
5-5700-IR	TechCut 5™ with Internal Reservoir, 100-240 V



Features:

- Unique fixturing design that allows use of either indexing or T-slot tables
- Soft-starting function for easing blade into cut/material
- Adjustable cut depth with high speed auto retraction
- Sample rotation for difficult, round or thick samples (requires rotation attachment)
- Touchpad switches to control all functions



The TechCut 5's unique design allows the use of either a T-slot table or X-axis indexing table, which are easily and quickly changed to support a wide range of fixture configurations.

The increased surface of the T-slot table is ideal for larger samples. The indexing tables provide quantitative positioning of the sample to the blade for precision sectioning applications.


This flexibility enables quick breakdown and setup of fixturing when dealing with a wide range of sample shapes and sizes.

Easy to Use



Featuring an intuitive control panel and backlit display, the TechCut 5™ is easy to operate, allowing greater efficiency among users. Select between imperial or metric units. All keypad buttons feature icon-based graphics that are easy to read, including Y-axis stage movement, RPM, force, coolant and sample rotation.

Automatic Sample Feed



During the cutting process, adaptive force control optimizes the feed rate through the sectioning process, allowing the operator to select force or feed as the preference for control.

As thicker or harder regions of the sample pass through the blade, the resistance will increase. The adaptive control lowers the feed rate as it passes through this region and raises it when it detects less resistance.

Cut depth can be set either with the pre-defined options, ranging from from 0.25" (6 mm) to 5" (125 mm), or by the user via the ADJ (adjustable) setting up to 4.25" (108 mm).

Jog functionality allows the operator to quickly move the sample into position near the blade once the sample is secured.

Manual Sample Feed



Manual Y-axis sample advancement using the hand wheel provides additional flexibility in positioning and feeding. The wheel is utilized to slowly feed delicate samples into the blade for minimal damage, allowing the user to "feel" the sectioning operation or to quickly position the sample prior to sectioning. It is also convenient for blade dressing.

Dual Arbor Flange




The TechCut 5™ accommodates blades between 3" (75 mm) and 8" (200 mm) in diameter, with either a 0.5" (12.7 mm) or 1.25" (32 mm) ID arbor hole.



When reversed, the outer flange features a step-arbor to allow use of blades with a 1.25" (32 mm) arbor hole.


Large Cutting Compartment



With a 5" (125 mm) clearance to the left-hand side of the blade and 6" (150 mm) to the right side, long and/or large parts are easily accommodated.

The T-slot table provides 2" (50 mm) of clearance beneath the spindle, and the X-axis tables provide 3" (75 mm) of clearance.

Blade Dressing




Loaded/clogged blades can easily be dressed to resume optimal cutting performance. The quick-slide vise or vertical clamp is used to hold the dressing stick, and then indexed manually into the blade using the hand wheel.

LED Illumination



LED illumination inside the cover allows easy observation of the sectioning process through the large viewing windows.

Wash-Down Hose



The wash-down hose is used to clean the cutting compartment. A quick-connect fitting on the end of the coolant nozzle bridge allows the hose to be connected for temporary use and disconnected for storage.

T-Slot Table



The T-slot table allows a wide range of clamping and fixturing options available to accommodate a variety of sample shapes and sizes.

Measuring 7.75" W x 8.25" D (195 x 210 mm), it is machined from aluminum billet and fitted with stainless tops for durability.

Item	Description
5-5330	T-Slot Table

Table Clamp



Flat samples up to 0.25" (6 mm) thick can easily be secured to the T-slot table with these height adjustable clamps.

Item	Description
5-5350	Table Clamp (Pk/4)


Vertical Adjusting Clamp



The vertical adjusting clamp can be positioned in any orientation throughout the T-slots in the table. The height adjustable lever allows various size samples up to 2.25" (57 mm) to be secured. Interchangeable pads allow handling of a variety of sample shapes. One (1) angled pad and one (1) flat pad are included.

Item	Description
5-4050	Vertical Adjusting Clamp
5-4050SPA	Swivel Pad Assortment Kit


Quick-Slide Vise



Made of hardened 17-4 stainless steel, the quick-slide vise features a lead-screw that disengages for quick and easy operation. Removable delrin inserts (included) are ideal for holding brittle or delicate samples. The vise mounts directly to the indexing or T-slot tables.

Item	Description
5-5525	Quick-Slide Stainless Steel Vise

X-Axis Indexing Tables




The X-axis indexing tables allow quantitative positioning for sectioning in precise locations. Two resolution options are offered: 25 μm (coarse) and 5 μm (fine). Coarse positioning is more applicable for metallurgical sectioning applications, where fine positioning is more suitable for tight geometry requirements usually found in microelectronic failure analysis and dicing applications.

Both feature 4" (100 mm) of travel and a 0.5" (12.7 mm) center-threaded mounting interface for maximum fixture positioning versatility. The easy-to-read indexing dial can be rotated to allow zeroing of the scale.

Item	Description
5-5320	X-Axis Indexing Table, 25 μm Resolution
5-5325	X-Axis Indexing Table, 5 μm Resolution


Saddle Clamps



These stainless steel saddle clamps provide single or dual clamping for cross-section cuts. They mount directly onto the indexing tables and can be spaced at different widths for versatility. V-notches on both sides of the clamp hold parts securely in place.

Item	Description
5-5555	Saddle Clamp (Pk/2)

Quick-Slide Vise




Made of hardened 17-4 stainless steel, the quick-slide vise features a lead-screw that disengages for quick and easy operation. Mounts directly to the indexing or T-slot tables.



Removable delrin inserts (included) are ideal for holding brittle or delicate samples. (Vise riser #5-5515 shown above.)

Item	Description
5-5525	Quick-Slide Stainless Steel Vise
5-5514	0.75" (19 mm) Vise Riser
5-5515	1.5" (38 mm) Vise Riser

Cam-Lock Fixture Adapter



The cam-lock fixture adapter offers a quick and convenient method for securing and detaching the #15-1045 multipurpose fixture to the indexing tables. The sectioned sample can then be transferred without disturbing orientation directly to the MultiPrep™ for grinding/polishing on the same plane.

Item	Description
5-5520	Cam-Lock Fixture Adapter
15-1045	Multipurpose Fixture

Rotational Cutting Attachment



The rotational cutting attachment, which mounts directly to the indexing or T-slot tables, enables sample rotation during a cut. Its benefits include:

- Increased sample thickness capacity
- Shortened cutting time
- Reduced blade/motor load
- Reduced heat on the sample

The sample can be rotated continuously or oscillated through the cut. A variety of fixtures are included (see page 10).

Item	Description
5-5745	Rotational Cutting Attachment with Holders

Dicing Rotation Stage



The dicing rotation stage (shown mounted on #5-5730 Z-axis stage and #5-5325 fine indexing table) allows 360° rotational positioning of flat substrates. Samples are mounted to the borosilicate glass insert using wax, which is secured with set-screws for quick attachment and removal. It can be mounted directly to the indexing tables, or to the Z-stage to control cutting depth.

Item	Description
5-5540	Dicing Rotation Stage, 360°/0.1° resolution

Teardrop Holder



The teardrop holder secures various sized mounts from 25 mm to 2" diameter. It is mounted to the indexing tables for precise location sectioning, and positions the mount at the center of the blade.

Item	Description
5-5560	Teardrop Holder for 25 mm - 2" mounts


Fixture Rotation Stage



The fixture rotation stage (shown with #5-5525 vise) allows angle cutting in 2° increments. A variety of fixtures can be mounted to provide rotational adjustment of the sample relative to the blade. It mounts directly to the indexing tables or Z-axis stage.

Item	Description
5-5070	Fixture Rotation Stage, 360°/2° resolution

Z-Axis Stage



The Z-axis stage (shown mounted on #5-5325 fine indexing table with #5-5540 dicing rotation stage) allows vertical positioning of the sample relative to the blade. It provides 2" (51 mm) of travel with 5 µm increments. The platform has a 0.5" center-threaded mounting interface for stage/fixture positioning versatility, and it mounts directly to the indexing tables.

Item	Description
5-5730	Z-Axis Stage, 5 µm increments

Multiposition Vise



The multiposition vise allows samples of up to 3" (76 mm) to be secured. It can be mounted directly to the indexing tables, or to the Z-axis or fixture rotation stages.

Item	Description
5-5510	Multiposition Vise

Coolant System Options

Internal Reservoir, 1.25 Gallon (4.7 L)

An internal 1.25-gallon (4.7 L) stainless steel reservoir sits beneath the cutting platform. An internal pump draws coolant from the reservoir through a filtered intake designed to prevent clogging.



Shown with optional #5-5580 Magnetic Particle Collector

External Reservoir, 7 Gallon (26.5 L)

An external, 7-gallon (26.5 L) stainless steel tank that sits below the TechCut 5™ provides added capacity for high-use applications where necessary maintenance of the coolant is less frequent. Power to the pump is supplied through the auxiliary outlet on the back panel. This system features a "catch-filter" for the drain that is easily accessible for cleaning (inset).




External Filtration System




With a reusable/cleanable 75 µm stainless steel filter, this filtration system is designed for applications where excessive debris is generated during sectioning operations (plastics, resins, composites, ferrous metals).

Item	Description
5-5590	External Filtration System

Universal T-Slot / X-Axis Table Accessories




5-5525
Quick-Slide Vise with delrin jaw inserts, 2" (51 mm) capacity




5-5745
Rotational Cutting Attachment

T-Slot Table Accessories




5-4050
Vertical Adjusting Clamp




5-5350
Table Clamp


X-Axis Table Accessories




5-5555
Dual Saddle Clamp set, with 0.25" (6 mm) risers




5-5510
Multiposition Vise, 3" (76 mm) capacity




5-5520
Cam-Lock Fixture Adapter




5-5560
Teardrop Holder, 25 mm - 2" mounts



5-5540
Dicing Rotation Stage




5-5730
Z-Axis Stage




5-5070
Fixture Rotation Stage


Rotational Attachment Fixtures




20-1501
Surface Mount
Fixture



20-1507
Mount Holder,
2" (51 mm)
diameter capacity



20-1515
Multipurpose
Fixture,
1" x 1.5"
(25 x 38 mm)
capacity



20-1502
Multipurpose
Fixture,
1.5" x 2"
(38 x 51 mm)
capacity


Cabinet/Stand




5-4230
Cabinet



Cabinet - Internal View with
External Coolant Reservoir
(5-5700-ER)




Cabinet (Front Elevation)



Cabinet (Side Elevation)

Dimensions



Technical Specifications

Blade Range	Diameter: 3" - 8" (76 - 200 mm) Arbor Hole: 0.5" or 1.25" (12.7 or 32 mm)
Flange Size	2" (51 mm) & 3" (76 mm) included
Sample/Cutting Capacity	6" W x 6" L x 2.5" H (152 x 152 x 64 mm)
Y-Axis Table	8" (203 mm) travel - manual or programmable cut-depth of 0.05" - 3" (1 - 76 mm) per minute
X-Axis Table	4" (102 mm) travel, 5 µm or 25 µm increments with zeroing dial
T-Slot Table	7.75" D x 8.25" W (195 x 210 mm) with four (4) channels
Force	Low / Medium / High
Motor Power	1.25 HP (950 W)
Blade RPM	500 - 5,000 (100 RPM increments)
Electronics	Microprocessor controlled
Operating Voltage	100-240 VAC, single phase 50/60 Hz
Display	LCD 4 x 16 - imperial or metric units
Cutting Action(s)	Linear or Linear + Rotation
Recirculating Coolant Capacity	Internal: 1.25 gallons (4.7 L) External: 7 gallons (26.5 L)
Safety	Electronic sensor/blade cut-off for cover, Emergency stop
CE	LVD Directive 2006/95/EC, EN 60204-1:2006+A1:2009, EMC Directive 2004/108/EC
Weight	147 lb. (67 kg)
Warranty	2 years
Includes	Operation manual, Country specific power cord

